

# Experiment 5 : Study Of Microwave Components

## Aim

To determine the frequency & wavelength in a rectangular waveguide working on TE<sub>10</sub> mode.

## Instruments/Equipments

1. Klystron Power Supply
2. Klystron tube with Klystron mounts
3. Isolator
4. Variable attenuator
5. Frequency meter
6. Slotted section
7. Tunable probe
8. oscilloscope
9. BNC cable

## Theory

Mode represents in wave guides as either TE<sub>mn</sub> / TM<sub>mn</sub>, Where TE-Transverse electric, TM-Transverse magnetic.  $m$  - Number of half wave length variation in broader direction.  $n$  - Number of half wave length variation in shorter direction.

$$\frac{\lambda_g}{2} = d_1 - d_2$$

Where  $d_1$  and  $d_2$  are the distance between two successive minima/maxima. It is having highest cut off frequency hence dominant mode. For dominant TE<sub>10</sub> mode in rectangular wave guide  $\lambda_0, \lambda_g, \lambda_c$  are related as below.

$$\frac{1}{\lambda_0^2} = \frac{1}{\lambda_g^2} + \frac{1}{\lambda_c^2}$$

Where  $\lambda_0$  is free space wave length,  $\lambda_g$  is guide wave length,  $\lambda_c$  is cutoff wave length. For TE<sub>10</sub> mode

$$\lambda_c = \frac{2a}{m}$$


Where  $m = 1$  in TE<sub>10</sub> mode and  $a$  is inner broad dimension of waveguide. The wavelength of the signal in an unbounded medium (air or vacuum), calculated as

$$\lambda_0 = c/f$$

Where  $c = 3 \times 10^8$  m/s is velocity of light and  $f$  is frequency. For propagation to occur, the operating free space wavelength must be less than the cutoff wavelength ( $\lambda_0 < \lambda_c$ )

## Procedure

1. Set up the components and equipments as shown in figure.
2. Set Mode selector switch to FM-Mode position with FM amplitude and FM frequency knob at mid position. Keep beam voltage control knob fully anticlockwise(minimum) and reflector voltage knob to fully clockwise(Maximum).
3. Fan should be kept in front of klystron
4. Switch on Fan



5. Switch On the klystron power supply and oscilloscope. Adjust the repeller voltage until a square wave on a DSO. Record the parameters beam voltage, beam current, repeller voltage correctly using the Mode Select switch on the Klystron Power Supply